

STERILE LINERS

- Easy to install
- Complete coverage of tray top and sides:
 - liners attach at tray bottom
- Liners protect tray and clearly identify **Clean** vs. **Dirty**
- Compatible with other transport trays



SafeCap saves us the cost of torn liners, positive ATP tests and additional cleaning.

GI Lab Manager
300 bed hospital
Southern California

SafeCap is easier to install and provides better patient and staff protection.

GI Director
Midwest Healthcare
System



 **CLEAN**



 **DIRTY**



TRANSPORT TRAYS AND LIDS

- Reusable, high temperature washable trays and lids
- Elegant design: no difficult to clean areas
- Ergonomic shelf for easy handling
- Solid, secure lid allows for stacking of trays
- Compatible with other transport carts

Easy to clean!



Meets SGNA[®], AORN[®], and AAMI[®] Guidelines

| Features | SafeCap | Competitive System | Bag System | Disposable Tray System |
|--|---------|--------------------|------------|------------------------|
| Reusable, high temperature washable trays and lids | ✓ | ✗ | ✗ | ✗ |
| Solid lid design allows for safe stacking of trays | ✓ | ✗ | ✗ | ✗ |
| Easy to read liner identification:  CLEAN  DIRTY | ✓ | ✓ | ✗ | ✗ |
| Liners cover trays sides and attach at the bottom of the tray | ✓ | ✗ | ✗ | ✗ |
| Custom sterile liners to ensure contamination kept in/out on top and sides | ✓ | ✗ | ✗ | ✗ |
| Easily view red and green liners on transport cart | ✓ | ✓ | ✗ | ✗ |
| Easy to clean and disinfect tray — no lip or crevices | ✓ | ✗ | ✗ | ✗ |
| Tray supports weight of endoscope | ✓ | ✓ | ✗ | ✗ |

| ORDERING INFORMATION | | |
|----------------------|--|---------------|
| Part Number | Product Description | Units Per Box |
| SC-1478 | Sterile Liner Set (clear, green, and red liner) | 100 |
| SC-CLR-100 | Clear Sterile Liners - 100 ea. individually packaged | 100 |
| SC-GRN-100 | Green Sterile Liners - 100 ea. individually packaged | 100 |
| SC-RED-100 | Red Sterile Liners - 100 ea. individually packaged | 100 |
| SC-152L | Clear Reusable Lid** | 1 |
| SC-249G | Gray High Temperature Tray* | 1 |
| SC-249L | Gray High Temperature Lid* | 1 |
| SC-CART | Six tray cart* 51"(h) x 27.5"(w) x 25.5" (d) | 1 |

*Maximum temp. 230°F/110°C

**Maximum temp. 180°F/82°C



SafeCap Transport Cart

- Easy access
- Clear visibility
- Locking wheels
- Contoured cart top
- Secure tray transport